

The Original Microbrush® Applicator

Consistent, reliable quality



4
Tip sizes

Ultrafine	0,5 mm
Superfine	1,0 mm
Fine	1,5 mm
Regular	2,0 mm

9
Colors



3
Packaging options

Dispenser Kit
Dispenser Refill
Tube



Item #	Description			
MPD	Microbrush® Plus Dispenser Kit	● yellow	fine	1,5 mm
Microbrush® Plus Applicator Dispenser Kit				
PR400KIT	Microbrush® Plus Applicator Dispenser Kit	assorted	regular	2,0 mm
PF400KIT	Microbrush® Plus Applicator Dispenser Kit	● pink ● yellow	fine	1,5 mm
PSF400KIT	Microbrush® Plus Applicator Dispenser Kit	● white	superfine	1,0 mm
Microbrush® Plus Applicator Refill				
PR400	Microbrush® Plus Applicator Refill	assorted	regular	2,0 mm
PR400BL	Microbrush® Plus Applicator Refill	● blue	regular	2,0 mm
PR400GR	Microbrush® Plus Applicator Refill	● green	regular	2,0 mm
PR400PE	Microbrush® Plus Applicator Refill	● pfirsich	regular	2,0 mm
PR400PU	Microbrush® Plus Applicator Refill	● violett	regular	2,0 mm
PF400	Microbrush® Plus Applicator Refill	● pink ● yellow	fine	1,5 mm
PF400PI	Microbrush® Plus Applicator Refill	● pink	fine	1,5 mm
PF400YE	Microbrush® Plus Applicator Refill	● yellow	fine	1,5 mm
PSF400	Microbrush® Plus Applicator Refill	● white	superfine	1,0 mm
PU400TE	Microbrush® Plus Applicator Refill	● teal	ultrafine	0,5 mm
PU400OR	Microbrush® Plus Applicator Refill	● orange	ultrafine	0,5 mm
Microbrush® Tube Applicators				
MRA400	Microbrush® Tube Applicators	assorted	regular	2,0 mm
MRB400	Microbrush® Tube Applicators	● blue	regular	2,0 mm
MRG400	Microbrush® Tube Applicators	● green	regular	2,0 mm
MRS400	Microbrush® Tube Applicators	● pfirsich	regular	2,0 mm
MRP400	Microbrush® Tube Applicators	● violett	regular	2,0 mm
MFA400	Microbrush® Tube Applicators	● pink ● yellow	fine	1,5 mm
MFP400	Microbrush® Tube Applicators	● pink	fine	1,5 mm
MFY400	Microbrush® Tube Applicators	● yellow	fine	1,5 mm
MSF400	Microbrush® Tube Applicators	● white	superfine	1,0 mm
MUT400	Microbrush® Tube Applicators	● teal	ultrafine	0,5 mm
MUO400	Microbrush® Tube Applicators	● orange	ultrafine	0,5 mm






Microbrush® X

- For even placement of material to inner surface of tooth wall
- Extended and extra thin for use in confined spaces
- Designed with stiffer fibers for scrubbing
- For endodontic treatment, also for fissures and bonding posts
- Patented bending portion is bendable up to 90 degrees



MPX
Microbrush® X Dispenser
Kit, black, 100 ct.

PX100
Microbrush® X Refill, black,
100 ct.







STIX64B
Micro-Stix™ Applicator,
64 ct., blue


Micro-Stix™ Applicator

- Ideal for veneers, crowns, brackets, inlays and onlays
- Applicator tip bendable
- Innovative adhesive formula
- Enables to securely hold restorations until placed securely in their proper place.
- Can be used multiple times during a single patient procedure




UltraBrush[®] Bristle Brush

- Bristle brush applicators for spreading and smoothing materials
- Bends close to tip for improved control
- “Stay-Put Bending” Technology - Maintain the applicator’s head at desired angle without rebound
- Chemical-resistant fibers designed to hold up to the harshest dental materials
- Available in two sizes: UltraBrush[®] regular ± 200 fibers, UltraBrush[®] fine ± 100 fibers
- Ideal for desensitizers, whitening solutions, fluoride, sealants and more



Item #	Description		
U2D	UltraBrush [®] Dispenser Kit	● green	regular
U2R200	UltraBrush [®] Refill	● yellow / green	regular
U1D	UltraBrush [®] Dispenser Kit	● blue	fine
U1R200	UltraBrush [®] Refill	● orange / blue	fine



dry tips®

Also available
DryTips for
ultimate moisture
control

Far outlasts cotton rolls and other absorbent pads,
and doesn't lose moisture when handled.

DryTips® advantage over traditional cotton rolls:

- Effective coverage of parotid orifice and buccal mucosa
- Polyethylen film for cutting protection
- Enables easier access to molar region
- Flexible, accomodates cheek movement
- Stays comfortably in place, yet is easily removed with water spray
- Reflective coating for increased visibility and to aid in retraction



NEW!
Reflective coating
now available.

IDEAL FOR

- Impressions
- Occlusal adjustments
- Orthodontics
- Fluoride applications
- Composites
- Cementation
- Fissure sealing
- Paint-on techniques

Item #	Description	Color	Size	Quantity
291542	DryTips®	white	large	50 pc.
291543	DryTips®	white	small	50 pc.
291726	Reflective DryTips®	● silver	large	50 pc.
291727	Reflective DryTips®	● silver	small	50 pc.

MB-GB-CAT/0724

We are happy to help!

info@microbrush.eu | +49 6221 43 45 442

YOUNG
INNOVATIONS